

Section 1. Product and Company Identification

Item Number: s216ba-1
 Common Name: Weigerts Hematoxylin Solution A
 Intended Use: In Vitro Diagnostic use. Laboratory Use Only
 IN CASE OF EMERGENCY, CONTACT: CHEMTREC (24HR) 800-424-9300

Manufacturer: Poly Scientific R&D Corp

70 Cleveland Ave
 Bay Shore, NY 11706
 631-586-0400
 polyrnd@polyrnd.com
 www.polyrnd.com

Section 2. Hazard Identification

225 Flammable Liquids Cat 2
 302 Acute toxicity, oral Cat 4
 315 Skin corrosion/irritation Cat 2
 319 Serious eye damage/eye irritation Cat 2A
 335 Specific target organ toxicity, single exposure; Respiratory tract irritation Cat 3
 370 Specific target organ toxicity, single exposure Cat 1



Danger

Highly flammable liquid and vapour. Harmful if swallowed. Causes serious eye irritation. May cause respiratory irritation. Causes damage to eyes, blood and CNS. Keep away from heat/sparks/open flames/hot surfaces. -No smoking. Do not breathe dust/fume/gas/mist/vapours/spray. Use only in well-ventilated area. In case of inadequate ventilation wear respiratory ventilation. Wash hands/skin thoroughly after handling. Do not eat, drink or smoke when using this product. Wear protective gloves/protective clothing/eye protection/face protection. Keep container tightly closed, locked up in well ventilated-area and cool. Ground/bond container and receiving equipment. Use explosion-proof electrical/ventilating/lighting/equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Take off contaminated clothing and wash before reuse. In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish. If skin irritation occurs: Get medical advice/attention. If eye irritation persists: get medical advice/attention. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. Rinse Mouth. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Wash skin with soap/water/shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER or doctor/physician if you feel unwell. If exposed: Call a POISON CENTER or doctor/physician. Dispose of contents/container to an approved waste disposal plant.

Hazards not otherwise covered by GHS: None

Section 3. Composition Information

Exposure Limits(A blank value indicates no information available)

Component	CAS#	PEL(mg/m ³)	STEL(mg/m ³)	CEIL(mg/m ³)	Concentration Range
Ethyl Alcohol	64-17-5		1,900.00		50-100%
isopropyl Alcohol	67-63-0	1,225.00	980.00		0-5%
Hematoxylin	517-28-2				0-5%
Methyl Alcohol	67-56-1	325.00	260.00		0-5%

Section 4. First Aid Measures

Eye Contact: Check for and remove contact lenses. Wash with large amounts of water for 15 minutes. Seek medical attention.
 Skin Contact: Remove contaminated clothing and shoes. Wash the affected area with large with soap and water. Seek medical attention
 Ingestion: Give two glasses of water to a conscious victim. Do not induce vomiting. Seek medical attention
 Inhalation: Move person to fresh air. If necessary give CPR; warning this could pose a risk of exposure to the rescue breather. Seek medical attention
 The most important known symptoms and effects are described in section 2 and/or section 11.

Section 5. Fire Fighting Measures

Extinguishing Media: Use Dry Chemical, Foam or Carbon Dioxide
 Special Fire and Explosion Remarks: N/A

Section 6. Accidental Release Measures

Spill Cleanup and Disposal Special: Eliminate Ignition sources. Take up spill with absorbent material
 Spill Cleanup: Take up spills with absorbant material and containerize for proper disposal. Use proper PPE as per section 8. Provide ventilation.

Section 7. Handling and Storage

Storage and Handling Special: N/A
 Storage and handling: Keep container tightly closed. Store in a cool, dry area and protect from physical damage

Section 8. Exposure Controls/Personal Protection

Personal Protective Equipment: Safety Glasses, Gloves, Vapor Respirator
 This information is provided as a guide but proper PPE can only be determined by the end user and their situation.
 Engineering Controls: Provide local exhaust ventilation to keep the airborne concentrations of vapors below their respective threshold limit values. Ensure that eyewash stations and safety showers are local to the work-station.

Section 9. Physical and Chemical Properties

Appearance	clear tan to brown	Evaporation Rate.....	N/A	Water Soluable?	Yes
Odor	Pleasant	Upper Flammability Limit (%):	N/A	Volatile Percent	100
Odor Threshold.....	N/A	Lower Flammability Limit (%):	N/A	Partition Coefficient.....	n-octanol/water: N/A
pH.....	N/A	Specific Gravity (@ 20C).....	.791	Auto Ignition Temp.	N/A
Melting Point.....	N/A	Vapor Pressure (mm Hg).....	52	Decomposition Temp.....	N/A
Boiling Point.....	N/A	Vapor Density (Air=1)	1.6	Viscosity	N/A
Flash Point (F) TCC	N/A	Relative Density.....	N/A		

Section 10. Stability and reactivity

Special Remarks on Stability: Stable

Special Remarks on Reactivity: N/A

Water Reactive: No

Section 11. Toxicological Information

Routes of Entry: Inhalation, Skin Absorption, Ingestion

Animal Toxicity: Ethyl Alc: Acute Oral (LD50): 7060 mg/kg (Rat); Acute Dermal(LD50): 500 mg/24hr (Rabbit); Methyl Alc: Acute Oral (LD50):5628 mg/kg (Rat); Acute Dermal (LD50) 500mg/24hr (Rabbit); Isopropyl Alc: Acute Oral (LD50):5045 mg/kg (Rat); Acute Dermal (LD50) 500mg (Rabbit); Hematoxylin: Acute Oral (TDLO) 400 gm/kg (Rat)

Human Toxic Effects: Target Organs: Respiratory system, skin, eyes, CNS, liver, blood and reproductive system

Potential Acute Health Effects: Hazardous in case of inhalation, eye contact, skin contact, ingestion

Potential Chronic Health Effects: Isopropyl Alcohol: IARC Code 3

Section 12. Ecological Information

Ecological Information: N/A

Section 13. Disposal Considerations

Waste Disposal: Dispose of in accordance with local, state and federal laws.

Section 14. Transport Information

DOT Identification: UN1170; Ethanol Solutions: 3; II

Section 15. Regulatory Information

State Regulations: New York Release reporting list: N/A

Sara Section 311 Reporting

Component	CAS#	Acute	Chronic	Fire	Pressure	Reactive	SARA302	SARA313	CERCLA	RCRA
Ethyl Alcohol	64-17-5	No	No	No	No	No	No	Yes	Yes	No
isopropyl Alcohol	67-63-0	No	No	No	No	No	No	Yes	No	No
Hematoxylin	517-28-2	No	No	No	No	No	No	No	No	No
Methyl Alcohol	67-56-1	No	No	No	No	No	No	Yes	Yes	No

Section 16. Other Information

Review Date: 12/17/2015

Reviewed by: ddi

MSDS Group Id: 70

Notice: This MSDS applies only to the material as packaged. If the material is altered by any means it may pose risks not mentioned here.

It is the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use.

While this MSDS is based on reliable technical data, Poly Scientific R&D Corp. assumes no responsibility for the completeness or accuracy of the information contained herein.

Section 1. Product and Company Identification

Item Number: s216bb-1
 Common Name: Weigerts Hematoxylin Solution B
 Intended Use: In Vitro Diagnostic use. Laboratory Use Only
 IN CASE OF EMERGENCY, CONTACT: CHEMTREC (24HR) 800-424-9300

Manufacturer: Poly Scientific R&D Corp

70 Cleveland Ave
 Bay Shore, NY 11706
 631-586-0400
 polyrnd@polyrnd.com
 www.polyrnd.com

Section 2. Hazard Identification

290 Corrosive to Metals Cat 1
 318 Serious eye damage/eye irritation Cat 1



Danger

May be corrosive to metals. Causes severe skin burns and eye damage.

Wash hands/skin thoroughly after handling. Wear protective gloves/protective clothing/eye protection/face protection. Immediately call a POISON CENTER or doctor/physician. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician. Absorb spillage to prevent material damage. Keep only in original container and tightly closed. Dispose of contents/container to an approved waste disposal plant

Hazards not otherwise covered by GHS: None

Section 3. Composition Information

Exposure Limits(A blank value indicates no information available)

Component	CAS#	PEL(mg/m3)	STEL(mg/m3)	CEIL(mg/m3)	Concentration Range
Hydrochloric Acid	7647-01-0	5.00		5.00	0-5%
Ferric Chloride	10025-77-1				0-5%

Section 4. First Aid Measures

Eye Contact: Check for and remove contact lenses. Wash with large amounts of water for 15 minutes. Seek medical attention.

Skin Contact: Remove contaminated clothing and shoes. Wash the affected area with large with soap and water. Seek medical attention

Ingestion: Give two glasses of water to a conscious victim. Do not induce vomiting. Seek medical attention

Inhalation: Move person to fresh air. If neccessary give CPR; warning this could pose a risk of exposure to the rescue breather. Seek medical attention

The most important known symptoms and effects are described in section 2 and/or section 11.

Section 5. Fire Fighting Measures

Extinguishing Media: Use Dry Chemical, Carbon Dioxide, Water Spray, Alcohol Foam

Special Fire and Explosion Remarks: N/A

Section 6. Accidental Release Measures

Spill Cleanup and Disposal Special: Warning! Corrosive! Wear protective clothing and respiratory equipment. Pick up with absorbent material.

Spill Cleanup: Take up spills with absorbant material and containerize for proper disposal. Use proper PPE as per section 8. Provide ventilation.

Section 7. Handling and Storage

Storage and Handling Special: N/A

Storage and handling: Keep container tightly closed. Store in a cool, dry area and protect from physical damage

Section 8. Exposure Controls/Personal Protection

Personal Protective Equipment: Safety Glasses, Gloves, Vapor Respirator

This information is provided as a guide but proper PPE can only be determined by the end user and their situation.

Engineering Controls: Provide local exhaust ventilation to keep the airborne concentrations of vapors below their respective threshold limit values. Ensure that eyewash stations and safety showers are local to the work-station.

Section 9. Physical and Chemical Properties

Appearance	clear yellow	Evaporation Rate.....	N/A	Water Soluable?	Yes
Odor.....	Strong	Upper Flammability Limit (%):	N/A	Volatile Percent.....	100
Odor Threshold.....	N/A	Lower Flammability Limit (%):	N/A	Partition Coefficient.....	n-octanol/water: N/A
pH.....	N/A	Specific Gravity (@ 20C).....	1.18	Auto Ignition Temp.....	N/A
Melting Point.....	N/A	Vapor Pressure (mm Hg).....	190	Decomposition Temp.....	N/A
Boiling Point.....	N/A	Vapor Density (Air=1).....	N/A	Viscosity	N/A
Flash Point (F) TCC.....	N/A	Relative Density.....	N/A		

Section 10. Stability and reactivity

Special Remarks on Stability: Stable

Special Remarks on Reactivity: N/A

Water Reactive: No

Section 11. Toxicological Information

Routes of Entry: Inhalation, Skin Absorption, Ingestion

Animal Toxicity: Hydrochloric Acid : Acute Oral (LD50): 900 mg/kg (Rat); Ferric Chloride: Acute Oral (LDLOO 900 mg/kg (Rat) I.V. (LDLO) 7 mg/kg (Rabbit)

Human Toxic Effects: Target Organs: Respiratory System, Eyes, Skin, GI Tract, Liver

Potential Acute Health Effects: Hazardous in case of inhalation, eye contact, skin contact, ingestion

Potential Chronic Health Effects: Hydrochloric Acid: IARC Code 3

Section 12. Ecological Information

Ecological Information: N/A

Section 13. Disposal Considerations

Waste Disposal: Dispose of in accordance with local, state and federal laws.

Section 14. Transport Information

DOT Identification: UN1789; Hydrochloric Acid: 8; II

Section 15. Regulatory Information

State Regulations: New York Release reporting list: Hydrochloric Acid, Ferric Chloride

Sara Section 311 Reporting

Component	CAS#	Acute	Chronic	Fire	Pressure	Reactive	SARA 302	SARA 313	CERCLA	RCRA
Hydrochloric Acid	7647-01-0	No	No	No	No	No	Yes	Yes	Yes	No
Ferric Chloride	10025-77-1	No	No	No	No	No	No	No	No	No

Section 16. Other Information

Review Date: 12/17/2015

Reviewed by: ddi

MSDS Group Id: 72

Notice: This MSDS applies only to the material as packaged. If the material is altered by any means it may pose risks not mentioned here.

It is the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use.

While this MSDS is based on reliable technical data, Poly Scientific R&D Corp. assumes no responsibility for the completeness or accuracy of the information contained herein.