

Section 1. Product and Company Identification

Item Number: s2816-1

Common Name: Acetic Acid 2% in Reagent Alcohol

Intended Use: In Vitro Diagnostic use. Laboratory Use Only

IN CASE OF EMERGENCY, CONTACT: CHEMTREC (24HR) 800-424-9300

Manufacturer: Poly Scientific R&D Corp

70 Cleveland Ave
Bay Shore, NY 11706
631-586-0400
polyrnd@polyrnd.com
www.polyrnd.com**Section 2. Hazard Identification**

225 Flammable Liquids Cat 2

302 Acute toxicity, oral Cat 4

314 Skin corrosion/irritation Cat 1A, B, C

370 Specific target organ toxicity, single exposure Cat 1



Danger

Highly flammable liquid and vapour. Harmful if swallowed. Causes severe skin burns and eye damage. Causes damage to eyes, blood and CNS.

Keep away from heat/sparks/open flames/hot surfaces.-No smoking. Do not breathe dust/fume/gas/mist/vapours/spray. Wash hands/skin thoroughly after handling. Do not eat, drink or smoke when using this product. Wear protective gloves/protective clothing/eye protection/face protection. Keep container tightly closed, locked up in well ventilated-area and cool. Ground/bond container and receiving equipment. Use explosion-proof electrical/ventilating/lighting/equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Take off contaminated clothing and wash before reuse. In case of fire: Use dry sand, dry chemical or alcohol- resistant foam to extinguish. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. Rinse mouth. Do not induce vomiting. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Wash contaminated clothing before reuse. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.. If exposed: Call a POISON CENTER or doctor/physician. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON CENTER or doctor/physician. Dispose of contents/container to an approved waste disposal plant.

Hazards not otherwise covered by GHS: None

Section 3. Composition Information

Exposure Limits(A blank value indicates no information available)

Component	CAS#	PEL(mg/m3)	STEL(mg/m3)	CEIL(mg/m3)	Concentration Range
ETHyl Alcohol	64-17-5		1,900.00		50-100 %
Isopropyl Alcohol	67-63-0	1,225.00	980.00		0-5%
Acetic Acid, Glacial	64-19-7	15.00	10.00		0-5%
Methyl Alcohol	67-56-1	325.00	260.00		0-5%

Section 4. First Aid Measures**Eye Contact:** Check for and remove contact lenses. Wash with large amounts of water for 15 minutes. Seek medical attention.**Skin Contact:** Remove contaminated clothing and shoes. Wash the affected area with large with soap and water. Seek medical attention**Ingestion:** Give two glasses of water to a conscious victim. Do not induce vomiting. Seek medical attention**Inhalation:** Move person to fresh air. If necessary give CPR; warning this could pose a risk of exposure to the rescue breather. Seek medical attention

The most important known symptoms and effects are described in section 2 and/or section 11.

Section 5. Fire Fighting Measures**Extinguishing Media:** Use Dry Chemical, Foam, or Carbon Dioxide**Special Fire and Explosion Remarks:** N/A**Section 6. Accidental Release Measures****Spill Cleanup and Disposal Special:** Warning! Flammable! Eliminate ignition sources. Take up spill with absorbent material.**Spill Cleanup:** Take up spills with absorbant material and containerize for proper disposal. Use proper PPE as per section 8. Provide ventilation.**Section 7. Handling and Storage****Storage and Handling Special:** N/A**Storage and handling:** Keep container tightly closed. Store in a cool, dry area and protect from physical damage**Section 8. Exposure Controls/Personal Protection****Personal Protective Equipment:** Safety Glasses, Gloves, Vapor Respirator

This information is provided as a guide but proper PPE can only be determined by the end user and their situation.

Engineering Controls: Provide local exhaust ventilation to keep the airborne concentrations of vapors below their respective threshold limit values. Ensure that eyewash stations and safety showers are local to the work-station.**Section 9. Physical and Chemical Properties**

Appearance	clear, colorless	Evaporation Rate	N/A	Water Soluable?.....	Yes
Odor	Slight vinegar	Upper Flammability Limit (%).....	N/A	Volatile Percent.....	100
Odor Threshold.....	N/A	Lower Flammability Limit (%).....	N/A	Partition Coefficient.....	n-octanol/water: N/A
pH.....	N/A	Specific Gravity (@20C).....	0.791	Auto Ignition Temp.....	N/A
Melting Point	N/A	Vapor Pressure (mm Hg).....	52	Decomposition Temp	N/A
Boiling Point.....	N/A	Vapor Density (Air=1)	1.6	Viscosity	N/A
Flash Point (F) TCC	N/A	Relative Density	N/A		

Section 10. Stability and reactivity**Special Remarks on Stability:** Stable**Special Remarks on Reactivity:** N/A**Water Reactive:** No

Section 11. Toxicological Information

Routes of Entry: Inhalation, Skin Absorption, Ingestion

Animal Toxicity: Ethyl Alcohol: Acute Oral (LD50) 7060 mg/kg (Rat) Acute Dermal 500mg/24hr (Rabbit); Methyl Alcohol: Acute Oral (LD50) 5628 mg/kg (Rat) Acute Dermal 500mg/24hr (Rabbit); Isopropyl Alcohol: Acute Oral (LD50) 5045 mg/kg (Rat) Acute Dermal 500mg/24hr (Rabbit); Glacial Acetic Acid: Acute Oral (LD50): 3310 mg/kg (Rat); Acute Dermal(LD50):

Human Toxic Effects: Target Organs: Eyes, Skin, Respiratory System, CNS, Liver, Blood, Reproductive System

Potential Acute Health Effects: Hazardous in case of inhalation, eye contact, skin contact, ingestion

Potential Chronic Health Effects: Isopropyl Alcohol: IARC Code 3

Section 12. Ecological Information

Ecological Information: N/A

Section 13. Disposal Considerations

Waste Disposal: Dispose of in accordance with local, state and federal laws.

Section 14. Transport Information

DOT Identification: UN1170; Ethanol Solutions; 3; II

Section 15. Regulatory Information

State Regulations: New York release reporting list: Acetic Acid

Component	CAS#	Sara Section 311 Reporting								CERCLA	RCRA
		Acute	Chronic	Fire	Pressure	Reactive	SARA302	SARA313			
Ethyl Alcohol	64-17-5	No	No	No	No	No	No	Yes	Yes	No	
Isopropyl Alcohol	67-63-0	No	No	No	No	No	No	Yes	Yes	No	
Acetic Acid, Glacial	64-19-7	No	No	No	No	No	No	No	Yes	No	
Methyl Alcohol	67-56-1	No	No	No	No	No	No	Yes	Yes	No	

Section 16. Other Information

Review Date: 6/14/2018

Reviewed by: tsc

MSDS Group Id: 93

Notice: This SDS applies only to the material as packaged. If the material is altered by any means it may pose risks not mentioned here.

It is the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use.

While this SDS is based on reliable technical data, Poly Scientific R&D Corp. assumes no responsibility for the completeness or accuracy of the information contained herein.