

Section 1. Product and Company Identification

Item Number: C826-1GL
 Common Name: Ortho Xylene
 Intended Use: In Vitro Diagnostic use. Laboratory Use Only
 IN CASE OF EMERGENCY, CONTACT: CHEMTREC (24HR) 800-424-9300

Manufacturer: Poly Scientific R&D Corp

70 Cleveland Ave
 Bay Shore, NY 11706
 631-586-0400
 polyrnd@polyrnd.com
 www.polyrnd.com

Section 2. Hazard Identification

226 Flammable Liquids Cat 3
 315 Skin corrosion/irritation Cat 2
 332 Acute toxicity, inhalation Cat 4
 361 Reproductive toxicity Cat 2



Warning

Flammable liquid and vapor.
 Causes skin irritation and serious eye irritation.
 May cause irritation to the respiratory system.
 May be fatal if swallowed and enters airways.

Obtain special instructions before use. Do not handle until all precautions have been read and understood. Keep away from heat/sparks/open flames/hot surfaces. No smoking. Use explosion-proof electrical/ventilating/lighting/equipment. Use only non-sparking tools. Take precautionary measures against static discharge. Wear protective gloves/protective clothing/eye protection/face protection. Avoid breathing dust/fume/gas/mist/vapours/spray. Use only in a well-ventilated area. Wash hands/skin thoroughly after handling. If exposed: Call a POISON CENTER or doctor/physician. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Call a POISON CENTER or doctor/physician if you feel unwell. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. Rinse mouth. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower. If skin irritation or rash occurs: Get medical advice/attention. Wash contaminated clothing before reuse. Store in well-ventilated place. Keep cool. Keep container tightly closed. Dispose of contents/container to an approved waste disposal plant.

Hazards not otherwise covered by GHS: None

Section 3. Composition Information

Exposure Limits (A blank value indicates no information available)

Component	CAS#	PEL(mg/m ³)	STEL(mg/m ³)	CEIL(mg/m ³)	Concentration Range
Xylene	95-47-6	150.00	100.00		50-100 %

Section 4. First Aid Measures

Eye Contact: Check for and remove contact lenses. Wash with large amounts of water for 15 minutes. Seek medical attention.
 Skin Contact: Remove contaminated clothing and shoes. Wash the affected area with large amount of soap and water. Seek medical attention.
 Ingestion: Give two glasses of water to a conscious victim. Do not induce vomiting. Seek medical attention.
 Inhalation: Move person to fresh air. If necessary give CPR; warning this could pose a risk of exposure to the rescue breather. Seek medical attention.
 The most important known symptoms and effects are described in section 2 and/or section 11.

Section 5. Fire Fighting Measures

Extinguishing Media: Dry Chemical, Carbon Dioxide, Water Spray, Alcohol Foam
 Special Fire and Explosion Remarks: N/A

Section 6. Accidental Release Measures

Spill Cleanup and Disposal Special: Eliminate all ignition sources. Pick up with absorbent vermiculite
 Spill Cleanup: Take up spills with absorbent material and containerize for proper disposal. Use proper PPE as per section 8. Provide ventilation.

Section 7. Handling and Storage

Storage and Handling Special: N/A
 Storage and handling: Keep container tightly closed. Store in a cool, dry area and protect from physical damage

Section 8. Exposure Controls/Personal Protection

Personal Protective Equipment: Splash Goggles, Gloves, Synthetic Apron, Vapor Respirator
 This information is provided as a guide but proper PPE can only be determined by the end user and their situation.
 Engineering Controls: Provide local exhaust ventilation to keep the airborne concentrations of vapors below their respective threshold limit values. Ensure that eyewash stations and safety showers are local to the work-station.

Section 9. Physical and Chemical Properties

Appearance	Evaporation Rate.....: N/A	Water Soluable?	Yes
Odor.....: Benzene like	Upper Flammability Limit (%): N/A	Volatile Percent	100
Odor Threshold.....: N/A	Lower Flammability Limit (%): N/A	Partition Coefficient.....:	n-octanol/water: N/A
pH.....: N/A	Specific Gravity (@ 20C).....: .87	Auto Ignition Temp.:	N/A
Melting Point.....: N/A	Vapor Pressure (mm Hg).....: 9	Decomposition Temp.....:	N/A
Boiling Point.....: N/A	Vapor Density (Air=1).....: 3.7	Viscosity	N/A
Flash Point (F) TCC.....: N/A	Relative Density.....: N/A		

Section 10. Stability and reactivity

Special Remarks on Stability: Stable

Special Remarks on Reactivity: N/A

Water Reactive: No

Section 11. Toxicological Information

Routes of Entry: Inhalation, Skin Absorption, Ingestion

Animal Toxicity: Acute Oral (LD50) 4300 mg/kg (Rat); Inhalation (LC50) 5000 ppm/4hr (Rat); Acute Dermal (LD50) 20000 mg/kg (Rabbit); Intraperitoneal (LD50) 1548 mg/kg (Mouse); Subcutaneous (LD50) 1700 mg/kg (Rat)

Human Toxic Effects: Target Organs: Eyes, Liver, Kidneys

Potential Acute Health Effects: Hazardous in case of inhalation, eye contact, skin contact, ingestion

Potential Chronic Health Effects: Xylene: IARC Code 3

Section 12. Ecological Information

Ecological Information: N/A

Section 13. Disposal Considerations

Waste Disposal: Dispose of in accordance with local, state and federal laws.

Section 14. Transport Information

DOT Identification: UN1307; Xylenes 3, III

Section 15. Regulatory Information

State Regulations: New York release reporting list: o-Xylene

Sara Section 311 Reporting

Component	CAS#	Acute	Chronic	Fire	Pressure	Reactive	SARA302	SARA313	CERCLA	RCRA
Xylene	95-47-6	No	No	No	No	No	No	Yes	Yes	No

Section 16. Other Information

Review Date: 12/17/2015

Reviewed by: ddi

MSDS Group Id: 304

Notice: This MSDS applies only to the material as packaged. If the material is altered by any means it may pose risks not mentioned here.

It is the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use.

While this MSDS is based on reliable technical data, Poly Scientific R&D Corp. assumes no responsibility for the completeness or accuracy of the information contained herein.