

## Section 1. Product and Company Identification

Item Number: C827-3LB  
 Common Name: Infiltrating/Embedding Paraffin Prills  
 Intended Use: In Vitro Diagnostic use. Laboratory Use Only  
 IN CASE OF EMERGENCY, CONTACT: CHEMTREC (24HR) 800-424-9300

Manufacturer: Poly Scientific R&D Corp

70 Cleveland Ave  
 Bay Shore, NY 11706  
 631-586-0400  
 polyrnd@polyrnd.com  
 www.polyrnd.com

## Section 2. Hazard Identification

Harmful If Swallowed or Inhaled

Avoid breathing dust or vapor. May be irritating to eyes, skin, and respiratory system. Wear safety goggles and rubber gloves to avoid contact. Wash thoroughly after handling. Keep container tightly closed.

EFFECTS OF EXPOSURE: Ingestion may cause nausea and abdominal pain. Mild irritant to skin and eyes.

TARGET ORGANS: None.

FIRST AID: Call a physician at once!

For Fire: Use extinguishing media appropriate for surrounding fire.

For Spill: Eliminate ignition sources. Pick up with absorbent material and containerize for proper disposal.

Hazards not otherwise covered by GHS: None

## Section 3. Composition Information

Exposure Limits(A blank value indicates no information available)

Component	CAS#	PEL(mg/m <sup>3</sup> )	STEL(mg/m <sup>3</sup> )	CEIL(mg/m <sup>3</sup> )	Concentration Range
Ammonium Hydroxide	1336-21-6	35.00	50.00		10-25%

## Section 4. First Aid Measures

Eye Contact: Check for and remove contact lenses. Wash with large amounts of water for 15 minutes. Seek medical attention.

Skin Contact: Remove contaminated clothing and shoes. Wash the affected area with large with soap and water. Seek medical attention

Ingestion: Give two glasses of water to a conscious victim. Do not induce vomiting. Seek medical attention

Inhalation: Move person to fresh air. If neccessary give CPR; warning this could pose a risk of exposure to the rescue breather. Seek medical attention

The most important known symptoms and effects are described in section 2 and/or section 11.

## Section 5. Fire Fighting Measures

Extinguishing Media: Dry Chemical, Carbon Dioxide, Water Spray, Alcohol Foam

Special Fire and Explosion Remarks: N/A

## Section 6. Accidental Release Measures

Spill Cleanup and Disposal Special: Pick up with absorbent vermiculite

Spill Cleanup: Take up spills with absorbant material and containerize for proper disposal. Use proper PPE as per section 8. Provide ventilation.

## Section 7. Handling and Storage

Storage and Handling Special: N/A

Storage and handling: Keep container tightly closed. Store in a cool, dry area and protect from physical damage

## Section 8. Exposure Controls/Personal Protection

Personal Protective Equipment: Safety Glasses, Gloves, Synthetic Apron

This information is provided as a guide but proper PPE can only be determined by the end user and their situation.

Engineering Controls: Provide local exhaust ventilation to keep the airborne concentrations of vapors below their respective threshold limit values. Ensure that eyewash stations and safety showers are local to the work-station.

## Section 9. Physical and Chemical Properties

Appearance .....	Evaporation Rate.....	Water Soluable? .....
Odor.....	Upper Flammability Limit (%):	Volatile Percent.....
Odor Threshold.....	Lower Flammability Limit (%):	Partition Coefficient.....
pH.....	Specific Gravity (@ 20C).....	Auto Ignition Temp.....
Melting Point.....	Vapor Pressure (mm Hg).....	Decomposition Temp.....
Boiling Point.....	Vapor Density (Air=1).....	Viscosity.....
Flash Point (F) TCC.....	Relative Density.....	

**Section 10. Stability and reactivity**

Special Remarks on Stability: Stable

Special Remarks on Reactivity: N/A

Water Reactive: No

**Section 11. Toxicological Information**

Routes of Entry: Inhalation, Skin Absorption, Ingestion

Animal Toxicity: Acute Oral (LD50) 350 mg/kg (Rat); Acute Oral (LD50) 43 mg/kg (Human); Inhalation (LCLO) 5000 PPM (Human); Intravenous (LDLO) 10 mg/kg (Rabbit)

Human Toxic Effects: Target Organs: Blood, Kidneys, Heart, CNS, Liver

Potential Acute Health Effects: Target Organs: Blood, Kidneys, Heart, CNS, Liver

Potential Chronic Health Effects: N/A

**Section 12. Ecological Information**

Ecological Information: N/A

**Section 13. Disposal Considerations**

Waste Disposal: Dispose of in accordance with local, state and federal laws.

**Section 14. Transport Information**

DOT Identification: Non Hazardous

**Section 15. Regulatory Information**

State Regulations: New York release reporting list: Ammonium Hydroxide

## Sara Section 311 Reporting

Component	CAS#	Acute	Chronic	Fire	Pressure	Reactive	SARA 302	SARA 313	CERCLA	RCRA
Ammonium Hydroxide	1336-21-6	No	No	No	No	No	No	No	Yes	No

**Section 16. Other Information**

Review Date: 12/17/2015

Reviewed by: ddi

MSDS Group Id: 305

Notice: This MSDS applies only to the material as packaged. If the material is altered by any means it may pose risks not mentioned here.

It is the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use.

While this MSDS is based on reliable technical data, Poly Scientific R&amp;D Corp. assumes no responsibility for the completeness or accuracy of the information contained herein.