

Section 1. Product and Company Identification

Item Number: s105-1
 Common Name: Acid Alcohol 3%
 Intended Use: In Vitro Diagnostic use. Laboratory Use Only
 IN CASE OF EMERGENCY, CONTACT: CHEMTREC (24HR) 800-424-9300

Manufacturer: Poly Scientific R&D Corp

70 Cleveland Ave
 Bay Shore, NY 11706
 631-586-0400
 polyrnd@polyrnd.com
 www.polyrnd.com

Section 2. Hazard Identification

225 Flammable Liquids Cat 2
 290 Corrosive to Metals Cat 1
 302 Acute toxicity, oral Cat 4
 315 Skin corrosion/irritation Cat 2
 318 Serious eye damage/eye irritation Cat 1
 370 Specific target organ toxicity, single exposure Cat 1



Danger

Highly flammable and liquid vapour. May be corrosive to metals. Harmful if swallowed. Causes skin irritation. Causes serious eye damage. Causes damage to organs. Keep away from heat/sparks/open flames/hot surfaces. -No smoking. Do not breathe dust/fume/gas/mist/vapours/spray. Wash hands/skin thoroughly after handling. Do not eat, drink or smoke when using this product. Wear protective gloves/protective clothing/eye protection/face protection. Keep container tightly closed, locked up in a dry, well ventilated-area, cool and in original container. Absorb spillage to prevent material damage. Use explosion-proof electrical/ventilating/lighting/equipment. Use only non-sparking tools. Take precautionary measures against static discharge. In case of fire: Use dry sand, dry chemical or alcohol-resistant foam to extinguish. IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell. Rinse mouth. Do not induce vomiting. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. Wash contaminated clothing before reuse. If skin irritation or rash occurs: Get medical advice/attention. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician. If exposed: Call a POISON CENTER or doctor/physician. Dispose of contents/container to an approved waste disposal plant.

Hazards not otherwise covered by GHS: None

Section 3. Composition Information

Exposure Limits(A blank value indicates no information available)

Component	CAS#	PEL(mg/m ³)	STEL(mg/m ³)	CEIL(mg/m ³)	Concentration Range
Isopropyl Alcohol	67-63-0	1,225.00	980.00		0-5%
Ethyl Alcohol	64-17-5		1,900.00		50-100 %
Hydrochloric Acid	7647-01-0	5.00		5.00	0-5%
Methyl Alcohol	67-56-1	325.00	260.00		0-5%

Section 4. First Aid Measures

Eye Contact: Check for and remove contact lenses. Wash with large amounts of water for 15 minutes. Seek medical attention.
 Skin Contact: Remove contaminated clothing and shoes. Wash the affected area with large with soap and water. Seek medical attention
 Ingestion: Give two glasses of water to a conscious victim. Do not induce vomiting. Seek medical attention
 Inhalation: Move person to fresh air. If necessary give CPR; warning this could pose a risk of exposure to the rescue breather. Seek medical attention
 The most important known symptoms and effects are described in section 2 and/or section 11.

Section 5. Fire Fighting Measures

Extinguishing Media: Dry Chemical Powder, Alcohol Foam, Carbon Dioxide
 Special Fire and Explosion Remarks: N/A

Section 6. Accidental Release Measures

Spill Cleanup and Disposal Special: Warning! Corrosive. Warning! Flammable! Eliminate ignitions. Take up spills with absorbent material.
 Spill Cleanup: Take up spills with absorbant material and containerize for proper disposal. Use proper PPE as per section 8. Provide ventilation.

Section 7. Handling and Storage

Storage and Handling Special: N/A
 Storage and handling: Keep container tightly closed. Store in a cool, dry area and protect from physical damage

Section 8. Exposure Controls/Personal Protection

Personal Protective Equipment: Gloves, Splash Goggles, Vapor Respirator, Apron
 This information is provided as a guide but proper PPE can only be determined by the end user and their situation.
 Engineering Controls: Provide local exhaust ventilation to keep the airborne concentrations of vapors below their respective threshold limit values. Ensure that eyewash stations and safety showers are local to the work-station.

Section 9. Physical and Chemical Properties

Appearance.....: clear colorless	Evaporation Rate.....: N/A	Water Soluable?.....: Yes
Odor.....: Strong	Upper Flammability Limit (%): N/A	Volatile Percent.....: 100
Odor Threshold.....: N/A	Lower Flammability Limit (%): N/A	Partition Coefficient.....: n-octanol/water: N/A
pH.....: N/A	Specific Gravity (@ 20C).....: .791	Auto Ignition Temp.....: N/A
Melting Point.....: N/A	Vapor Pressure (mm Hg).....: 52	Decomposition Temp.....: N/A
Boiling Point.....: N/A	Vapor Density (Air=1).....: 1.6	Viscosity.....: N/A
Flash Point (F) TCC.....: N/A	Relative Density.....: N/A	

Section 10. Stability and reactivity

Special Remarks on Stability: Stable

Special Remarks on Reactivity: N/A

Water Reactive: No

Section 11. Toxicological Information

Routes of Entry: Inhalation, Skin absorption, Ingestion

Animal Toxicity: Ethyl Ether: Acute Oral (LD50): 9750 mg/kg (Rat); Acute Dermal (LD50): 360 mg/kg (Rabbit); Acetone: : Acute Oral (LD50): 9.75 mg/kg (Rat); Acute Dermal(LD50): 20 mg/kg (Rabbit); Hydrochloric Acid: Acute Oral (LD50) 900 mg/kg (Rat) Intraperitoneal (LD50) 1449 mg/kg (Mouse)

Human Toxic Effects: Target Organs: Respiratory system, skin, eyes, CNS, liver, blood and reproductive system

Potential Acute Health Effects: Hazardous in case of eye, skin contact, inhalation, ingestion

Potential Chronic Health Effects: Isopropyl Alcohol: IARC Code 3; Hydrochloric Acid: IARC Code 3

Section 12. Ecological Information

Ecological Information: N/A

Section 13. Disposal Considerations

Waste Disposal: Dispose of in accordance with local, state and federal laws.

Section 14. Transport Information

DOT Identification: UN1170; Ethanol Solutions: 3; II

Section 15. Regulatory Information

State Regulations: New York release reporting list: Hydrochloric Acid

Sara Section 311 Reporting

Component	CAS#	Acute	Chronic	Fire	Pressure	Reactive	SARA302	SARA313	CERCLA	RCRA
Isopropyl Alcohol	67-63-0	No	No	No	No	No	No	Yes	No	No
Ethyl Alcohol	64-17-5	No	No	No	No	No	No	Yes	Yes	No
Hydrochloric Acid	7647-01-0	No	No	No	No	No	Yes	Yes	Yes	No
Methyl Alcohol	67-56-1	No	No	No	No	No	No	Yes	Yes	No

Section 16. Other Information

Review Date: 12/17/2015

Reviewed by: ddi

MSDS Group Id: 5

Notice: This MSDS applies only to the material as packaged. If the material is altered by any means it may pose risks not mentioned here.

It is the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use.

While this MSDS is based on reliable technical data, Poly Scientific R&D Corp. assumes no responsibility for the completeness or accuracy of the information contained herein.