

Section 1. Product and Company Identification

Item Number.: s2153-1
 Common Name.: Poly Mount (Xylene Based)
 Intended Use : In Vitro Diagnostic use. Laboratory Use Only
 IN CASE OF EMERGENCY, CONTACT: CHEMTREC (24HR) 800-424-9300

Manufacturer.: Poly Scientific R&D Corp.

70 Cleveland Ave
 Bay Shore NY 11706

polyrnd@polyrnd.com

Section 2. Hazard Identification

226 Flammable Liquids Cat 3
 304 Aspiration hazard Cat 1
 312 Acute toxicity,dermal Cat 4
 315 Skin corrosion/irritation Cat 2
 320 Serious eye damage/eye irritation Cat 2B
 332 Acute toxicity,inhalation Cat 4
 335 Specific target organ toxicity, single exposure; Respiratory tract irritation Cat 3
 350 Carcinogenicity Cat 1A, 1B
 361 Reproductive toxicity Cat 2
 370 Specific target organ toxicity, single exposure Cat 1

**Warning**

Flammable liquid and vapour. May be fatal if swallowed and enters airways. Harmful if inhaled and in contact with skin. Causes skin, eye and respiratory irritation. May cause cancer. Suspected of damaging fertility or the unborn child. Causes damage to CNS, Kidneys and Liver.

Obtain special instructions before use. Do not breathe dust/fume/gas/mist/vapours/spray. Wash hands/skin thoroughly after handling. Keep away from heat/sparks/open flames/hot surfaces.-No smoking. Avoid breathing dust/fume/gas/mist/vapours/spray. Wash skin thoroughly after handling. Do not eat, drink or smoke when using this product. Contaminated work clothing should not be allowed out of the workplace. Use only in well-ventilated area. In case of inadequate ventilation wear respiratory ventilation. Wear protective gloves/protective clothing/eye protection/face protection. Keep container locked up, tightly closed up in well ventilated-area and cool. Use explosion proof electrical/ventilating/lighting/equipment. Take precautionary measures against static discharge. Take off contaminated clothing and wash before reuse. In case of fire: Use dry sand, dry chemical or alcohol- resistant foam to extinguish. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with plenty of soap water/shower. Wash contaminated clothing before reuse. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician. If exposed: Call a POISON CENTER or doctor/physician. If skin irritation or rash occurs: Get medical advice/attention. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Dispose of contents/container to an approved waste disposal plant.

Hazards not otherwise covered by GHS: None

Section 3. Composition Information

Exposure Limits(A blank value indicates no information available) The specific chemical identity and/or exact percentage (concentration) of composition has been withheld as a trade secret.

Component	CAS#	PEL(mg/m3)	STEL(mg/m3)	CEIL(mg/m3)	Concentration Range
Tricresyl Phosphate	1330-78-5		0.10		<5%
Xylene	1330-20-7	651.00	434.00		>50%

Section 4. First Aid Measures

Eye Contact : Check for and remove contact lenses. Wash with large amounts of water for 15 minutes. Seek medical attention.
 Skin Contact : Remove contaminated clothing and shoes. Wash the affected area with large with soap and water. Seek medical attention
 Ingestion : Give two glasses of water to a conscious victim. Do not induce vomiting. Seek medical attention
 Inhalation : Move person to fresh air. If necessary give CPR; warning this could pose a risk of exposure to the rescue breather. Seek medical attention
 The most important known symptoms and effects are described in section 2 and/or section 11.

Section 5. Fire Fighting Measures

Extinguishing Media.: Dry Chemical, Carbon Dioxide, Water Spray, Alcohol Foam
 Special Fire and Explosion Remarks ...: N/A

Section 6. Accidental Release Measures

Spill Cleanup and Disposal Special ...: Danger! Flammable! Eliminate all ignition sources. Pick up with absorbent vermiculite
 Spill Cleanup.: Take up spills with absorbant material and containerize for proper disposal. Use proper PPE as per section 8. Provide ventilation.

Section 7. Handling and Storage

Storage and Handling Special...: N/A
 Storage and handling.: Keep container tightly closed. Store in a cool, dry area and protect from physical damage

Section 8. Exposure Controls/Personal Protection

Personal Protective Equipment ...: Splash Goggles, Gloves, Synthetic Apron, Vapor Respirator
 This information is provided as a guide but proper PPE can only be determined by the end user and their situation.
 Engineering Controls.: Provide local exhaust ventilation to keep the airborne concentrations of vapors below their respective threshold limit values. Ensure that eyewash stations and safety showers are local to the work-station.

Section 9. Physical and Chemical Properties

Appearance.....: Colorless liquid	Evaporation Rate.....: N/A	Water Soluable?.....: Yes
Odor.....: Benzene like	Upper Flammability Limit (%): N/A	Volatile Percent.....: 100
Odor Threshold.....: N/A	Lower Flammability Limit (%): N/A	Partition Coefficient.....: n-octanol/water: N/A
pH.....: N/A	Specific Gravity (@20C).....: .87	Auto Ignition Temp.....: N/A
Melting Point.....: N/A	Vapor Pressure (mm Hg).....: 9	Decomposition Temp.....: N/A
Boiling Point.....: N/A	Vapor Density (Air=1).....: 3.7	Viscosity.....: N/A
Flash Point (F) TCC.....: N/A	Relative Density.....: N/A	

Section 10. Stability and reactivity

Special Remarks on Stability.: Stable

Special Remarks on Reactivity.: N/A

Water Reactive.: No

Section 11. Toxicological Information

Routes of Entry.: Inhalation, Skin Absorption, Ingestion

Animal Toxicity.: Xylene: Acute Oral (LD50) 4300 mg/kg (Rabbit) Acute Dermal (LD50) 2000 mg/kg (Rabbit); Ethyl Benzene: N/A; Tricresyl Phosphate: Acute Oral (LD50) 3900 mg/kg (Mouse) Acute Dermal (LD50) 1500 mg/kg (Cat); Elvacite: N/A

Human Toxic Effects.: Target Organs: Eyes, skin, respiratory system, central nervous system, gastrointestinal tract, blood

Potential Acute Health Effects.: Hazardous in case of inhalation, eye contact, skin contact, ingestion

Potential Chronic Health Effects.: Xylene: IARC Code 3; Ethyl Benzene: IARC Code 2B

Section 12. Ecological Information

Ecological Information.: N/A

Section 13. Disposal Considerations

Waste Disposal.: Dispose of in accordance with local, state and federal laws.

Section 14. Transport Information

DOT Identification.: UN1307, Xylenes, 3, III

Section 15. Regulatory Information

State Regulations.: New York release reporting list: Xylene (Mixed)

Component	CAS#	Sara Section 311 Reporting						SARA302	SARA313	CERCLA	RCRA
		Acute	Chronic	Fire	Pressure	Reactive					
Tricresyl Phosphate	1330-78-5	No	No	No	No	No	No	No	No	No	
Xylene	1330-20-7	No	No	No	No	No	Yes	Yes	Yes	No	

Section 16. Other Information

Review Date : 3/15/2023

Reviewed by : Admin

MSDS Group Id.: 197

Notice: This SDS applies only to the material as packaged. If the material is altered by any means it may pose risks not mentioned here.

It is the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use.

While this SDS is based on reliable technical data, Poly Scientific R&D Corp. assumes no responsibility for the completeness or accuracy of the information contained herein.