

## Section 1. Product and Company Identification

Item Number: s2609-1  
 Common Name: Glutaraldehyde 2.5% In PBS pH 7.4  
 Intended Use: In Vitro Diagnostic use. Laboratory Use Only  
 IN CASE OF EMERGENCY, CONTACT: CHEMTREC (24HR) 800-424-9300

Manufacturer: Poly Scientific R&D Corp

70 Cleveland Ave  
 Bay Shore, NY 11706  
 631-586-0400  
 polyrnd@polyrnd.com  
 www.polyrnd.com

## Section 2. Hazard Identification

301 Acute toxicity, oral Cat 3  
 314 Skin corrosion/irritation Cat 1A, B, C  
 317 Sensitizer, Skin Cat 1  
 331 Acute toxicity, inhalation Cat 3  
 334 Sensitizer, respiratory Cat 1



Danger

Toxic if swallowed or if inhaled.  
 Causes severe skin burns and eye damage.  
 May cause an allergic skin reaction.  
 May cause allergy or asthma symptoms or breathing difficulties if inhaled.

Do not breathe dust/ fume/ gas/ mist/ vapours/ spray. Wash hands/skin thoroughly after handling. Do not eat, drink or smoke when using this product. Use only in a well-ventilated area. Contaminated work clothing should not be allowed out of the workplace. Wear protective gloves/ protective clothing/ eye protection/ face protection. In case of inadequate ventilation wear respiratory protection. IF SWALLOWED: Immediately call a POISON CENTER. Rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/ shower. If skin irritation or rash occurs: Get medical advice/ attention. Wash contaminated clothing before reuse. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. Immediately call a POISON CENTER or doctor/physician. If experiencing respiratory symptoms: Call a POISON CENTER or doctor/physician. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Store in a well-ventilated place. Keep container tightly closed and locked up. Dispose of contents/ container to an approved waste disposal plant.

Hazards not otherwise covered by GHS: None

## Section 3. Composition Information

Exposure Limits (A blank value indicates no information available)

Component	CAS#	PEL(mg/m <sup>3</sup> )	STEL(mg/m <sup>3</sup> )	CEIL(mg/m <sup>3</sup> )	Concentration Range
Glutaraldehyde	111-30-8		0.25		0-5%

## Section 4. First Aid Measures

Eye Contact: Check for and remove contact lenses. Wash with large amounts of water for 15 minutes. Seek medical attention.  
 Skin Contact: Remove contaminated clothing and shoes. Wash the affected area with large with soap and water. Seek medical attention  
 Ingestion: Give two glasses of water to a conscious victim. Do not induce vomiting. Seek medical attention  
 Inhalation: Move person to fresh air. If necessary give CPR: warning this could pose a risk of exposure to the rescue breather. Seek medical attention  
 The most important known symptoms and effects are described in section 2 and/or section 11.

## Section 5. Fire Fighting Measures

Extinguishing Media : Dry Chemical, Water Spray, Alcoholic Foam, Carbon Dioxide  
 Special Fire and Explosion Remarks : N/A

## Section 6. Accidental Release Measures

Spill Cleanup and Disposal Special : Pick up with absorbent vermiculite  
 Spill Cleanup : Take up spills with absorbant material and containerize for proper disposal. Use proper PPE as per section 8. Provide ventilation.

## Section 7. Handling and Storage

Storage and Handling Special : Refrigerate  
 Storage and handling: Keep container tightly closed. Store in a cool, dry area and protect from physical damage

## Section 8. Exposure Controls/Personal Protection

Personal Protective Equipment : Face Shield, Gloves, Synthetic Apron  
 This information is provided as a guide but proper PPE can only be determined by the end user and their situation.  
 Engineering Controls: Provide local exhaust ventilation to keep the airborne concentrations of vapors below their respective threshold limit values. Ensure that eyewash stations and safety showers are local to the work-station.

**Section 9. Physical and Chemical Properties**

Appearance .....: clear colorless	Evaporation Rate .....: N/A	Water Soluable? .....: Yes
Odor .....: Sharp	Upper Flammability Limit (%) :: N/A	Volatile Percent .....: 100
Odor Threshold .....: N/A	Lower Flammability Limit (%) :: N/A	Partition Coefficient .....: n-octanol/water: N/A
pH.....: N/A	Specific Gravity (@ 20C) .....: 1.131	Auto Ignition Temp. ....: N/A
Melting Point .....: N/A	Vapor Pressure (mm Hg) .....: 15	Decomposition Temp .....: N/A
Boiling Point .....: N/A	Vapor Density (Air=1) .....: N/A	Viscosity .....: N/A
Flash Point (F) TCC .....: N/A	Relative Density .....: N/A	

**Section 10. Stability and reactivity**

Special Remarks on Stability : Stable

Special Remarks on Reactivity : N/A

Water Reactive: No

**Section 11. Toxicological Information**

Routes of Entry: Inhalation, Skin Absorption, Ingestion

Animal Toxicity: Acute Oral (LD50) 134 mg/kg (Rat); Acute Dermal (LD50) 2560 mg/kg (Rabbit); Eye 1mg (severe irritation) (Rabbit); Skin 6mg/3 days (severe irritation) (Human)

Human Toxic Effects: Target Organs: Eyes, skin, respiratory system

Potential Acute Health Effects : Hazardous in case of inhalation, eye contact, skin contact, ingestion

Potential Chronic Health Effects : Respiratory

**Section 12. Ecological Information**

Ecological Information : N/A

**Section 13. Disposal Considerations**

Waste Disposal : Dispose of in accordance with local, state and federal laws.

**Section 14. Transport Information**

DOT Identification : Non Hazardous

**Section 15. Regulatory Information**

State Regulations: New York Release reporting list: N/A

Sara Section 311 Reporting

Component	CAS#	Acute	Chronic	Fire	Pressure	Reactive	SARA 302	SARA 313	CERCLA	RCRA
Gluteraldehyde	111-30-8	No	No	No	No	No	No	No	No	No

**Section 16. Other Information**

Review Date: 12/17/2015

Reviewed by: ddi

MSDS Group Id : 61

Notice: This MSDS applies only to the material as packaged. If the material is altered by any means it may pose risks not mentioned here.

It is the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use.

While this MSDS is based on reliable technical data, Poly Scientific R&D Corp. assumes no responsibility for the completeness or accuracy of the information contained herein.