

## Section 1. Product and Company Identification

Item Number: S294-1  
 Common Name: Zenker Fluid Modified Contains Zinc Chloride  
 Intended Use: In Vitro Diagnostic use. Laboratory Use Only  
 IN CASE OF EMERGENCY, CONTACT: CHEMTREC (24HR) 800-424-9300

Manufacturer: Poly Scientific R&D Corp

70 Cleveland Ave  
 Bay Shore, NY 11706  
 631-586-0400  
 polyrnd@polyrnd.com  
 www.polyrnd.com

## Section 2. Hazard Identification

272 Oxidising liquids; Oxidising solids Cat 2  
 301 Acute toxicity, oral Cat 3  
 314 Skin corrosion/irritation Cat 1A, B, C  
 330 Acute toxicity, inhalation Cat 1, 2  
 340 Germ cell mutagenicity Cat 1A, 1B  
 350 Carcinogenicity Cat 1A, 1B  
 360 Reproductive toxicity Cat 1A, 1B  
 370 Specific target organ toxicity, single exposure Cat 1



Danger

May intensify fire; Oxidizer. Toxic if swallowed. Causes severe skin burns and eye damage. Fatal if inhaled. May cause genetic defects. May cause cancer. May damage fertility or the the unborn child. Causes damage to organs.

Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. Keep away from heat/sparks/open flames/hot surfaces/clothing/combustible materials. -No smoking. Do not breathe dust/fume/gas/mist/vapours/spray. Wash hands/skin thoroughly after handling. Do not eat, drink or smoke when using this product. Wear protective gloves/protective clothing/eye protection/face protection. Do not get in eyes or on skin/clothing. In case of inadequate ventilation wear respiratory protection. Keep container tightly closed, locked up in well ventilated-area and cool. Keep/store away from clothing/combustible materials. Take precaution to avoid mixing with combustibles. In case of fire: Use water to extinguish. Do not use dry chemicals or foams for extinction. CO2 or Halon may provide limited control. IF SWALLOWED: Immediately call a POISON CENTER or doctor/physician. Rinse mouth. Do not induce vomiting. IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Wash skin with soap/water/shower. Wash contaminated clothing before reuse. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER or doctor/physician. If exposed: Call a POISON CENTER or doctor/physician. IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing. If having respiratory symptoms: Immediately call a POISON CENTER or doctor/physician. Dispose of contents/container to an approved waste disposal plant.

Hazards not otherwise covered by GHS: None

## Section 3. Composition Information

Exposure Limits(A blank value indicates no information available)

Component	CAS#	PEL(mg/m3)	STEL(mg/m3)	CEIL(mg/m3)	Concentration Range
Zinc Chloride	7646-85-7	2.00	1.00		0-5%
Potassium Dichromate	7780-50-9		0.05		0-5%

## Section 4. First Aid Measures

Eye Contact: Check for and remove contact lenses. Wash with large amounts of water for 15 minutes. Seek medical attention.

Skin Contact: Remove contaminated clothing and shoes. Wash the affected area with large with soap and water. Seek medical attention

Ingestion: Give two glasses of water to a conscious victim. Do not induce vomiting. Seek medical attention

Inhalation: Move person to fresh air. If necessary give CPR; warning this could pose a risk of exposure to the rescue breather. Seek medical attention

The most important known symptoms and effects are described in section 2 and/or section 11.

## Section 5. Fire Fighting Measures

Extinguishing Media: Dry Chemical, Carbon Dioxide, Water Spray, or Foam

Special Fire and Explosion Remarks: N/A

## Section 6. Accidental Release Measures

Spill Cleanup and Disposal Special: Pick up spill with absorbent vermiculite

Spill Cleanup: Take up spills with absorbant material and containerize for proper disposal. Use proper PPE as per section 8. Provide ventilation.

## Section 7. Handling and Storage

Storage and Handling Special: N/A

Storage and handling: Keep container tightly closed. Store in a cool, dry area and protect from physical damage

## Section 8. Exposure Controls/Personal Protection

Personal Protective Equipment: Safety Glasses, Gloves, Vapor Respirator

This information is provided as a guide but proper PPE can only be determined by the end user and their situation.

Engineering Controls: Provide local exhaust ventilation to keep the airborne concentrations of vapors below their respective threshold limit values. Ensure that eyewash stations and safety showers are local to the work-station.

## Section 9. Physical and Chemical Properties

Appearance .....	red orange	Evaporation Rate.....	N/A	Water Soluable? .....	Yes
Odor.....	N/A	Upper Flammability Limit (%):	N/A	Volatile Percent .....	N/A
Odor Threshold.....	N/A	Lower Flammability Limit (%):	N/A	Partition Coefficient.....	n-octanol/water: N/A
pH.....	N/A	Specific Gravity (@ 20C).....	N/A	Auto Ignition Temp. ....	N/A
Melting Point.....	N/A	Vapor Pressure (mm Hg).....	N/A	Decomposition Temp.....	N/A
Boiling Point.....	N/A	Vapor Density (Air=1) .....	N/A	Viscosity .....	N/A
Flash Point (F) TCC.....	N/A	Relative Density.....	N/A		

## Section 10. Stability and reactivity

Special Remarks on Stability: Stable

Special Remarks on Reactivity: N/A

Water Reactive: No

## Section 11. Toxicological Information

Routes of Entry: Inhalation, Skin Absorption, Ingestion

Animal Toxicity: Potassium Dichromate: Acute Oral (LD50) 190 mg/kg (Mouse); Sodium Sulfate: Acute Oral (LD50) 5989 mg/kg (Mouse); Zinc Chloride: Acute Oral (LD50) 350 mg/kg (Rat)

Human Toxic Effects: Target Organs: Respiratory System, Eyes, Liver, Kidneys, Skin, GI System, GI Tracts

Potential Acute Health Effects: Hazardous in case of inhalation, eye contact, skin contact, ingestion

Potential Chronic Health Effects: Potassium Dichromate: NTP Code 1

## Section 12. Ecological Information

Ecological Information: N/A

## Section 13. Disposal Considerations

Waste Disposal: Dispose of in accordance with local, state and federal laws.

## Section 14. Transport Information

DOT Identification: UN1840; Zinc Chloride, Solution; 8; III

## Section 15. Regulatory Information

State Regulations: New York release reporting list: Zinc Chloride

## Sara Section 311 Reporting

Component	CAS#	Acute	Chronic	Fire	Pressure	Reactive	SARA302	SARA313	CERCLA	RCRA
Zinc Chloride	7646-85-7	No	No	No	No	No	No	Yes	Yes	No
Potassium Dichromate	7780-50-9	No	No	No	No	No	No	Yes	Yes	No

## Section 16. Other Information

Review Date: 12/17/2015

Reviewed by: ddi

MSDS Group Id: 121

Notice: This MSDS applies only to the material as packaged. If the material is altered by any means it may pose risks not mentioned here.

It is the user's responsibility to develop proper methods of handling and personal protection based on the actual conditions of use.

While this MSDS is based on reliable technical data, Poly Scientific R&amp;D Corp. assumes no responsibility for the completeness or accuracy of the information contained herein.